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Fabrikant

BAAT Medical Products B.V.
Havenkade 300
7555 RJ Hengelo, The Netherlands

This instruction is valid only for the 'OTNI IntraMedullary Femur Prosthesis (OFF)', an implantable system for transcutaneous and intramedullary fixation of an exoprosthesis in patients with an amputation.

Warning

Before using any medical device distributed by OTNI Implants BV, the surgeon shall carefully pay attention to the following recommendations, as well as specific information about the device being used and that is displayed in the product sheet technical information and / or the surgical technique. Lack or failure to observe the information contained within the present instruction for use exonerates the Manufacturer from any liability whatsoever.

Description of the system

OTNI Ossointegration Femur Prosthesis (OFF) is an intramedullary, transfemoral and transcutaneous bone fixation system of the lower limb. It is an implantable modular system that can be used to directly connect an exoprosthesis to the lower limb. It is composed of the femoral, stump or above-knee amputees by means of the OTNI HELI adapter or the OPL GV adapter or the OTNI 17r L Modular Ossointegrationconnector.

OTNI Ossointegration Femur Prosthesis (OFF) is a dual cone bone fixation system comprising of the following components:

- Cone Femoral Stem equipped with OTNI Proximal Screw, an intramedullary stem allows the direct connection by means of ossointegration into the femoral bone
- OTNI Healing Plug, a temporary healing distal plug prevents of tissue ingrowth into the stem
- OTNI DC Adapter, a transcutaneous double cone adapter connects the intramedullary stem to an exoprosthesis
- OTNI Locking Screw, an internal locking screw for the connection of the transcutaneous double cone adapter to the stem

Intended Use

The OTNI OFF is intended to provide a fixed support for the connection of an artificial limb prosthesis to the residual femoral bone after transfemoral amputation, in all those cases of complications after, or contraindications to the use of conventional socket connections. The OTNI Ossointegrated Femur Prosthesis provides a bone anchorage point, to support the bone in cases of treatment of traumatic injuries or in corrective orthopedic surgery, as well as to be replaced by an exoprosthesis because of an above-knee amputation. The OFF is a non-active surgical implant designed for long-term implantation inside the human body. It is intended to be used by (orthopaedic-) surgeons with good knowledge of the specific operative technique, in a standard orthopaedic environment, for skeletally mature patients.

Indications for use

- Trans-femoral amputation.
- Complications after, or contraindications to the use of conventional socket prostheses.

Contraindications

- Inflammatory or septic, acute or chronic or local or distal systemic processes.
- Insufficient bone quality in order to seek appropriate and adequate press-fit anchoring and adequate ossointegration of intramedullary stem (eg severe osteoporosis, osteopenia).
- Vascular, musculoskeletal and neurological disorders.
- Bone metabolism disorders.
- Long-term corticosteroid or chemotherapy treatments.
- Patients with psychiatric disorders or mental instability or patients who are unwilling or unable to follow the instructions for rehabilitation and aftercare indicated by the doctor.
- Smoag, alcohol abuse, drug use.
- BMI > 25
- Pregnancy, lactation.
- Skeletally immature.

General information and precautions for safe use of the system

Before the OTNI OFF system is used, it is always recommended that all details of the surgical technique and the instructions for use are understood by the surgeon. The OTNI OFF system has been validated (usability testing) when implanted with a specific instrumentation set.

- The Manufacturer validates a specific instrumentation set only for the fixation and optional removal of the OTNI OFF system. The Manufacturer disclaims any responsibility for the consequences arising from the use of tools and / or on the patient's body by third parties.
- A careful inspection of the device before its implantation must be carried out; there must be no damage present (scratches, spots, notches or deformations) or any dirt or signs of use. Also, the specific surgical instruments that are subject to wear and tear can lead to complications upon implantation. The Manufacturer and the Distributor recommends careful monitoring before using the specific instrumentation set, as well as for any of the OTNI OFF components with a visible damage or a damaged or damaged, do not use it. OTNI OFF components with a porous titanium coating must be used with care to prevent damage to the coating itself. Thus, the coated parts must not come into contact with materials other than the original packaging, surgical gloves.

- It is necessary to verify whether the primary and secondary sterile packaging of OTNI OFF components are intact and the expiration date of the package is necessary.
- The size of OTNI OFF system (intramedullary stem and transcutaneous double cone adapter) is selected by the surgeon based on preoperative planning based on standard radiographs and CT. The surgical procedure must be done with available complete range of all types and sizes of the system components and with the specific instruments set.

- The OTNI OFF intramedullary stem must be implanted with sufficient press-fit to ensure primary stability.
- During the press-fit implantation of the stem careful attention must be paid to prevent the occurrence of cracks or fractures.
- The use of the healing plug is optional, and the indication for its application is according to surgeon's discretion.
- The transcutaneous stoma of the stump should be placed so that there is enough space between the skin and the transcutaneous double cone adapter, in such way subcutaneous stoma moisture can easily drain out, and facilitate a good epithelialization of the stoma wall surrounding the transcutaneous adapter.
- For the placement of transcutaneous double cone adapter, the healing plug must be removed and all surfaces must be clean and dry.

- The length (size) of the transcutaneous double cone adapter must be that one for which the skin is in correspondence with the cylindrical part of the adapter and the distal more taper is completely outside with respect to the stump stoma (see Surgical Technique).
- The length of the transcutaneous double cone adapter is dependent on the soft tissues and is such that there is at least 2 cm of TiNn coating penetrates through the skin
- The tightening torque of the internal locking screw is 10 Nm.
- Three months after full load of the OTNI OFF system check the torque of the internal locking screw.
- The surgeon has full discretion to determine when one or more components of the OTNI OFF system must be removed or replaced.

- The TiNn coating of some components acts as an insulating layer for preventing release of metal ions from the underlying metal. The time that this insulating layer works is unknown and cannot be guaranteed, the risks associated with the release of the metal ions and the wear of the components cannot be excluded.
- The OTNI OFF components are single use and MUST NEVER BE REUTILIZED. It is not allowed to reuse an implantable component, that has previously been implanted in the body of the patient or a third person, or has been in contact with body fluid or tissue of a third person.
- It is advisable to observe any and all additional information present in the Surgical Technique and on the components packaging labels. Complications or other problems derived from improper indications or surgical techniques, choices for use of incompatible components, improper use of instruments are exclusively of surgeon's responsibility and not of the manufacturer.
- The patient is responsible for stoma cleaning at least twice daily.

Possible side effects

The possible consequences of the implantation of an OTNI OFF system are:

- General risks associated with anesthesia;
- General risks associated with major surgery;
- Cardiovascular complications including deep vein thrombosis, thrombophlebitis, pulmonary embolism, wound hematoma and avascular necrosis;
- Nerve damage, neuromas;
- Intra-osseous lesion, including fissures, perforations, peri-prosthetic fractures or stem removal;
- Abnormal stump pain (acute, chronic, psychosomatic);
- Heterotopic ossification;
- Bacterial colonization of the stoma or soft tissue area with or without signs or symptoms of infection;
- Complications in the stoma area including edema, swelling, tissue irritation including delay of the stoma epithelialization;
- Superficial infections of the stoma and / or soft tissue;
- Deep intramedullary infections, osteitis, osteomyelitis, septic arthritis;
- Muscle contractures, enthesitis, tendinitis;
- Lossing, fracture, or mechanical failure of prosthetic components due to excessive load (excessive physical activity, sport, inattention or contraindications) or non-physiological stress (fall and / or traumatic accidents) or due to incorrect choice of components (eg, improper size of stem and transcutaneous adapter);
- Falling in case of breakage of the prosthetic components;
- Possible excessive tension on the soft tissues of the stump as skin, scars or muscles;
- Peri-prosthetic bone resorption and bone loss, osteonecrosis;
- Retrograde stem migration;
- Lack or insufficient ossointegration of the intramedullary stem, with subsequent loosening, due to a g.poor bone quality, inadequate primary stability (press fit), micro-movements, early mobilization and loading;
- Allergic reaction to the implant materials or the emission of metal ions by the implant;
- Complications due to modular prosthetic connections (metallosis, adverse biological reactions, osteolysis due to particles and release of metal ions due to corrosion, abrasion and wear component interface)

- Revision surgery to replace one or more parts before the end of its useful life;
- Revision surgery for a permanent explant of a OTNI OFF;
- Revision surgery with amputation of the stump at a higher (more proximal) level;
- Detection by security surveillance systems based on electromagnetic fields and interference with clinical analysis or instrumental analyses based on magnetic fields;

Pre-operative planning

Inaccurate pre-surgical planning can lead to disappointing results, or a system failure. Pre-operative planning requires a proper selection of patients on the basis of indication and contra-indication criteria. The surgeon must evaluate the stump with respect to the soft tissue quality and the ipsilateral hip functions. The surgeon should also define the most appropriate size of the intramedullary stem to be implanted based on length, diameter and morphology of the residual femur on plain radiological images.

Information for the patient

The surgeon must provide the patient with all the information regarding the surgical procedure, limitations of the implant, postoperative treatment and aftercare. It is responsibility of the surgeon to provide information about the postoperative treatment, with concern to the field of rehabilitation treatment, clinical follow-up, information on cleaning the OTNI OFF system and soft tissue hygiene.

- The surgeon must give detailed information to the patient about potential general risks associated with surgical procedures, possible side effects and eventual intrinsic limitations of the system, as well as on preventive measures that should be taken to reduce any risk.
- It is important that the surgeon informs the patient about the following:
 - Requirement to the patient of a full commitment in attending rehabilitation and a certain amount of clinical visits;
 - The operation and maintenance of the OTNI OFF system
 - Amount of loads and movements allowed during postoperative rehabilitation.

The patient must have understood all the information prior to the implantation of the OTNI OFF system and must be motivated to undergo the procedure. The patient must strictly follow all the conditions established by the surgeon concerning the rehabilitation (modality and duration) frequency of clinical follow-up, cleaning the OTNI OFF system, the cleaning of the stoma and must inform the surgeon of unforeseen consequences and / or complications.

It is advisable the surgeon informs the patient about the following instructions for use, cleaning, hygiene and maintenance of the OTNI OFF system:

- Twice a day clean the skin and soft tissue around the stoma with clean water and mild soap, and then dry in the air.

Sterility

Each part of the OTNI OFF system is packed separately and delivered in a sterile package. It is advisable to check on the product label of the component whether the part is sterile.

The manufacturer provides the implantable components in the 'Sterile' package which should be maintained protected until the moment of implantation. Before use check always:

- The expiration date of the sterility (year / month) on the product label;
- Visually check whether the secondary packaging and labeling are intact.
- Visually check that the sterile primary packaging is intact and there are no cracks, holes or other damage signs present on it. The implantable devices are sterilized by irradiation with a dose of 25 kGy.

Materials used in implantable devices,

The label for each unit displays the data, and the type of materials / coating that has been used. Manufacturing materials are:

- Ti6Al7Nb ISO5832-11for intramedullary femoral stem, dual cone transcutaneous adapter and locking screw
- Ti6Al4V ISO5832-3 for proximal screw of the intramedullary femoral stem
- Ti CP ISO5832-2 for the coating of the intramedullary femoral stem
- UHMWPE ISO5834-2for healing plug and plug of the internal self locking screw
- TiNn, for coating of the parts of intramedullary femoral stem, proximal screw, and dual cone transcutaneous adapter

Lifetime

The OTNI OFF system is designed to anchor the exoprosthesis directly into the remaining femur in patients with transfemoral amputation, wherein the mechanical forces at full load can be transmitted for a limited time. In general, the functional duration of an exoprosthesis bone fixation system is limited because it is exposed to unavoidable wear and tear, aging and complications that may lead to surgical re-intervention with the possible removal or replacement of the system.

Storage conditions

This device shall be stored under conditions specified to prevent deterioration by sunlight, moisture, temperature or other conditions. Storage conditions shall be monitored and recorded periodically, where appropriate.

MRI safety

A patient with this implant/device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm (a higher value for the spatial gradient magnetic field may apply if properly calculated).
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (per pulse sequence).

In all cases, the Health Care Professional is responsible for MR Conditions, MR Imaging quality and patient safety. Any safety issues or major image artefacts should be reported.

MRI-Related Heating

In non-clinical testing, comparable devices produced a temperature rise of less than or equal to 6.0 degrees C using an MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15-minutes (per pulse sequence) of scanning in a 3-Tesla MR system.

Artifact

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant/device. In some cases, the artifact size relative to the size of the implant or device may be indicated.

Attention: Contact the manufacturer of this implant/device for further information, as needed

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